


# KING AIR 90

Spectrum Aeromed is a leading designer and provider of customized air medical and ambulance equipment.

The company works closely with air medical and transport crews to equip all types of fixed wing and rotor wing aircraft with high-specialized and innovative life support solutions.




## King Air 90 Setup Overview:

- Quick Change capability can be installed or removed from the aircraft in less than 15 minutes and requires minimal modifications to the aircraft.
- High quality, durable design with cost-effective installation.
- Built aircraft and patient care specific.
- Ability to secure a variety of medical equipment mounts, multiple IV mounts, accessory tables and patient loading accessories.
- Ability to secure and transport infants with the Incubator Transport Deck.
- Ability to house storage cabinets.



## System Features

- Fully self-contained and lightweight
- The Module System allows up to 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter.
- AC outlets, DC outlet and dual pneumatic outlets.
- Custom stretchers are specifically designed for patient comfort.
- Aircraft is improved for installation of 2800 and 20 Series Module Systems.
- The module systems may be Electrical or Non-Electrical based on customer preference.



1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)