


# PILATUS PC-12


Spectrum Aeromed is recognized as a world leader in creating customized and durable, highly functional air medical transport equipment that complies with aviation regulations. All engineering and design of our air ambulance life support systems comes from working hand-in-hand with air medical personnel worldwide.



*Single 2800 Series System Install with Cabinets*



*Dual 2800 Series System Install*



*Triple 2800 Series System Install*

## System Overview:

- Quick change capability can be installed or removed from aircraft in approx. 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

## System Features & Options:

- Fully self-contained and lightweight
  - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
  - AC outlets, DC outlet and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV Mounts



Single 2800 Series System with Med Wall and Cabinets



Dual 2800 Series System Install



Stretcher and Base Install

1815 23rd Avenue North • Fargo, ND 58102

701.235.5478 • Toll Free: 1.800.753.4340 • Fax: 701.235.5486

Email: [Aeromed@Spectrum-Aeromed.com](mailto:Aeromed@Spectrum-Aeromed.com) • [www.Spectrum-Aeromed.com](http://www.Spectrum-Aeromed.com)