

# HAWKER 1000 INSTALLATION

## DUAL PATIENT LIFE SUPPORT SYSTEM

### System Overview:

- A more permanent configuration to specific needs and specific auxiliary medical equipment
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories



*Dual Dedicated System Install*

### System Features:

- Dedicated base and overhead panel
  - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
  - AC outlets, DC outlet and dual pneumatic outlets
  - Multiple storage drawers
- Custom stretcher specifically designed for patient comfort
- Multiple IV and medical equipment mounting options on the overhead panels
- 4 drawer storage cabinet



*Dedicated Overhead Panels*