

# INFANT TRANSPORT DECK

## AIRBORNE 20H INCUBATOR

### Transport System:

- Self contained critical care system for infants
  - The deck has provisions to mount an incubator on one end and a cylinder housing for oxygen and air on the opposite end.
- Quick-Change Capability
  - Can be installed or removed from aircraft in less than 10 minutes and requires minimal modifications to aircraft



*System Installed in a Cessna Excel*

### System Benefits Include:

- High quality system
- Light weight and durable design
- Cost-effective installation
- The ITS Deck without oxygen cylinders weights only 43 lbs



*System Installed in a Bell 407*