

# JETSTREAM 3201 INSTALLATION

## MODULAR LIFE SUPPORT SYSTEM

### System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

### System Features:

- Fully self-contained and lightweight
  - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
  - AC outlets, DC outlet and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV Mounts

### Spectrum Aeromed System Options:

- [2800 Series System](#)
- [20 Series System](#)

### Potential Installations:


*\*Contact a Sales Representative for further Install Details*




*Single 2800 Series System Install with Stretcher Bridge*



*Dual 2800 Series System Install*



*Dual 2800 Series System Install*



*Dual 2800 Series System Install with Medical Cabinets*