

KING AIR 200 INSTALLATION

DUAL PATIENT LIFE SUPPORT SYSTEM

System Overview:

- A more permanent configuration to specific needs and specific auxiliary medical equipment
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

System Features:

- Dedicated base and overhead
 - Two 3500 liter oxygen tanks, three air compressors, dual vacuum pumps and two 1000 watt inverter
 - AC outlets and pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV and medical equipment mounting options on the overhead panel
- Three storage drawers in the base
- Two recessed areas for suction canisters

Potential Installations:


**Contact a Sales Representative for further Install Details*



Dual Dedicated System Install



Dedicated Overhead Panels



Dual Dedicated System with Dedicated Overhead Panels