

LEAR 35 INSTALLATION

MODULAR LIFE SUPPORT SYSTEM

System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

System Features:


- Fully self-contained and lightweight
 - 3500 liters of oxygen, dual air compressors, vacuum pump and 1000 watt inverter
 - AC outlets, DC outlet and dual pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV Mounts

Spectrum Aeromed System Options:


- [2800 Series System](#)
- [20 Series System](#)
- [2200 Series System](#)

Potential Installations:

**Contact a Sales Representative for further Install Details*



Single 20 Series System Install




Single 20 Series System Install without Overhead




Single 20 Series System Install



Single 2200 Series System Install



Single 2800 Series System Install



Single 2800 Series System Install