

System Overview:

- Quick change capability can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to the aircraft
- High quality, durable design and cost-effective installation
- Built aircraft-specific and customer-specific
- Allows for a variety of custom medical mounts and patient loading accessories

System Features:


- Fully self-contained and lightweight
 - 3500 liters of oxygen, air compressor, vacuum pump and 500 watt inverter
 - AC outlets and pneumatic outlets
- Custom stretcher specifically designed for patient comfort
- Multiple IV Mounts
- 3 drawer storage cabinet

Spectrum Aeromed System Options:

- [2800 Series System](#)

Potential Installations:

**Contact a Sales Representative for further Install Details*



Single 2800 Series System Install



Medical Cabinet

Single 2800 Series System Install