

# Short Box Module 2000

Designed for Aircraft with Limited Cabin Length



## Spectrum Aeromed's Short Box Module 2000

This system was especially designed for Lear 24 operators who want to leave the potty seat installed. The stretcher and top can be positioned to cantilever over the front and/or aft end of the base. Shown with a custom stretcher bridge which secures a variety of medical equipment/monitors in an accessible place for crew.

### Short Box Base Features:

- Fully self-contained, lightweight.
- Can be installed or removed from the aircraft in 10 minutes.
- Base features 3500 liters of oxygen.
- 1000-Watt inverter.
- Dual air compressors and vacuum pump; and electrical outlets.

**The Overhead** is moveable along the entire length of the base and contains:

- Set of pneumatic outlets.
- Exam light that is moveable along the entire length of base.