



2800 Series Life Support System with Patient Loader Attached

Benefits include:

- High Quality, Light Weight, Durable, Cost-Effective Installation.
- Quick Change Capability-System can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to aircraft.
- Fully Self Contained Critical Care System with 3500 L oxygen, dual air compressors, vacuum pump, inverter, 115-volt AC electrical outlets, 28 volt DC outlet, dual pneumatic outlets, exam light.
- Fully Adjustable Stretcher with expandable armrests and adjustable headrest, 2" thick foam pad with bacteriostatic vinyl covering.
- FAA Certified and installed under a STC.