

2800 Series Life Support System

Installed in a Pilatus PCXII



Benefits Include:

- **Patients are easily loaded or unloaded with the Patient Loading System.**
- **Quick Change Capability** can be installed or removed from aircraft in less than 15 minutes and requires minimal modifications to aircraft.
- **Fully Self Contained Critical Care System** with 3500 L oxygen, dual air compressors, vacuum pump, inverter, 115-volt AC electrical outlets, 28 volt DC outlet, dual pneumatic outlets, exam light.
- **Fully Adjustable Stretcher** with expandable arm rests and adjustable head rest, 2" thick foam pad with bacteriostatic vinyl covering.
- **High Quality, Light Weight, Durable Design and Cost-Effective Installation.**